ATECO RI code: 72.19.09- Research and Experimental Development in different Natural Sciences and Engineering's fields.

DETERMINATION OF THE DECONTAMINATION PROCESS OF
SURFACES BY THE APPLICATION OF THE HIGH LEVEL
DISINFECTANT *ADANTIUM PLUS* WITH *HYGIENIO*
MACHINE

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COMPLIANCE STATEMENT

In relation to the conducted study, it is stated that the following report describes in faithfully and truthfully way the procedures we adopted and the results we obtained.

Tools and procedures are certified in accordance with:

- A) ACCREDIA- UNI EN ISO/ IEC 17025:2005 " General requirements for the testing and calibration Laboratories competence" (Attachment III).
- B) IQNet / DQS- Certification ISO 9001:2008 (Attachment IV).
- C) HEALTHCARE INSTITUTE-Approval of the production and/or sterilization quality assurance system (Attachment V).
- D) UNI EN ISO 4833-1:2013 : Horizontal Method for the microorganisms count. Part 1: Colonies count at 30 °C with the deep insemination technique.
- E) ISO/TS 19036:2006/Amd 1:2009 Measurement uncertainty for low counts.
- F) UNI EN ISO 15883-1 e 2:2009- Washing and disinfection equipment.

The study was coordinated by the working group:

U. Sciamannini, Ph.D	_____	Date:	29/10/2014
	(<i>Study Director</i>)		

A. Sciamannini, A.D.	_____	Date:	29/10/2014
	(<i>General Administrator</i>)		

In collaboration with:

<i>Dr. S. Saveri</i>	<i>Senior Chemist FIS&DM Srl</i>
----------------------	--------------------------------------

<i>Dr L. Alzalamira</i>	<i>Biologist Jacaroni Diagnostic Centre</i>
-------------------------	---

<i>Dr. L.Molinari</i>	<i>Biologist Jacaroni Diagnostic Centre</i>
-----------------------	---

1. CLIENT DATA.....	4
2. TEST LABORATORY.....	4
3. TEST RELATED DATE.....	4
4. REFERENCE DOCUMENTS.....	5
4.1 International law.....	5
5. PURPOSE.....	5
6 TEST DATA.....	5
6.1 HYGIENIO device description.....	5
6.2 ADANTIUM PLUS disinfectant description.....	6
6.3 Sampling plane.....	6
7. MATERIALS AND METHODS USED.....	7
7.1 Culture fields and tools of bacterial determination.....	7
8. OPERATIVE PROCEDURES FOR THE MICROBIOLOGICAL COUNT.....	9
8.1 Execution methods of the biological analysis.....	9
8.2 Test execution mode.....	9
8.3 Colonies count mode.....	9
9. RESULTS.....	11
10. CONCLUSIONS.....	12
11. ATTACHMENTS.....	13

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3. TEST RELATED DATE

Object of the study is the determination of the disinfection effectiveness of the *ADANTIUM PLUS* disinfectant through the use of the *HYGIENIO sanitizing*.

Responsible laboratories of the execution of the tests are ACCREDIA certificate (Accreditation number: 0820 Rev.2 Deadline: 2016-05-05- Attachment III).

Tests were conducted on the 14th October 2014 at the production office of FIS & DM Srl.

The tests of bacterial contamination has been performed on flat surfaces of Inox steel.

4. REFERENCE DOCUMENTS

4.1 *International law.*

For the protocol described in this relation, reference was made to the international law below specified.

- **Farmacopea Europea** – Current edition
- **UNI EN ISO 4833:2004** – Horizontal method for microorganisms count – Colonies count at 30 °C technique
- **UNI EN 15883-1:2009** – Washing and disinfection tools - Part 1: General requirements, terms, definitions and tests
- **UNI EN 15883-2:2009** - Washing and disinfection tools - Parte 2: Requirements and tests for washing and disinfection tools for chirurgic tools, anaesthetic equipment, hollow bodies, floors, containers, utensils, glassware, etc., using thermal disinfection.

5. PURPOSE

The validation consists in the valuation of the efficiency of the washing/decontamination process, obtained by the synergic use of the “ADANTIUM PLUS” disinfectant, vaporized through the use of HYGIENIO disinfectant device in order to locate and estimate the presence/absence of vital microorganisms in certain condition of culture .

6 TEST DATA

6.1 *HYGIENIO device description*

Hygienio is an innovative system for the hot application of the disinfectant, it quickly allows disinfecting surfaces and material of every kind at a high level and with the reuse oh spaces and equipment.

Thanks to the HYGIENIO system it is possible to obtain a two-phase mixture (disinfectant + steam) that in contact with a generic surface condenses, allowing a continuous and homogeneous disinfectant treatment. The steam jets, full of disinfectant particles, invade, sticking to them, floors, walls, furniture and objects of every shape and material and no-solution of continuity material, that is, without leaving micro-spaces where bacteria, fungus and mildew that would normally continue to proliferate.

6.2 ADANTIUM PLUS disinfectant description

ADANTIUM PLUS is a modern, concentrated synthetic preparation for simultaneous high-level cleansing and disinfection aimed at completely decontaminating Biomedical Devices in hospitals, according to current safety laws.

Adantium® Plus contains a new synthetic salt which, alone or resuspended with formula, is able to effectively hinder the nucleasic activity of the RNase present in plasma, preserving the added synthetic RNA from degradation (see the study of Prof. Pistello on the clinical validation).

Adantium® Plus, furthermore protects the cellular RNA from degradation through RNase A. In conclusion, on the basis of the experimental proof obtained and considering the high resistance to the physical-chemical inactivation of the RNase, Adantium® Plus has an effective inhibitory role on these nucleases which, in a biological environment, has an important role in curbing/inactivating viral replication.

Active ingredients: Diphedac®, Complex surfactant (consisting of Alkyl-dimethyl-benzyl ammonium chloride.), Solvents Isopropylic Alchol and Surfactants.

6.3 Sampling plane

To conduct the bacteriological test, have been effected:

- A sample through sterile swab, after contamination of the surfaces;
- A sample through sterile swab, successive at HYGIENIO/Adantium Plus treatment, after contact with the contaminated surfaces;

7. MATERIALS AND METHODS USED.

7.1 Culture fields and tools of bacterial determination

Culture fields used for the detection of the contaminant charge:

Field

Type of field: Plate Count Agar (PCA)
Manufacturing company: Biogenetics
Site: 118480/073
Deadline: 01/12/15

atcc Shackles used for the contamination

Salmonella ATCC 14028
Escherichia coli ATCC 25922
Staphylococcus aureus ATCC 25923
Listeria monocytogenes ATCC 11

Every field already cited is sterile. The preparation of the fields has been conducted following the producer indications and/or according to the referring method.

Features of the equipment used for the test:

Incubator at 30°C

Inventory number:01
builder: TERMOSTABIL Kt
Model: K2 TE
Measure and test field: 20°C-70°C, used at 30°C ± 1°C

Fridge - freezer

Inventory number:04
Builder: INDESIT
Modell: TAAN 3V
Measurement and test field: fridge used from +2°C / +8°C

Fridge

Inventory number:06
Builder: KELVINATOR
Model: KG 101
Measurement and test field: fridge used from +2°C / +5°C

Colonies count

Inventory number:16
Builder: STUART
Model: SC6
Measurement and test field: na

Logger date (inside the incubator)

Inventory number:09
Serial number:0526-0116
Builder: Escort Data Logging Systems Ltd, New Zealand
Model: 62 D 32
Measurement and test field: used at 30°C ± 1°C

Logger date (inside the fridge 4)

Inventory number:21
Serial number: MI-BH-358-082
Builder: Escort Data Logging Systems Limited
Model: MI-0E-D-2-L
Measurement and test field: - 40°C/+65°C used between 2°C e 5°C

8. OPERATIVE PROCEDURES FOR THE MICROBIOLOGICAL COUNT

8.1 Execution methods of the biological analysis

A 600 x 400 mm inox steel surface has been contaminated with the ATCC bacterial inoculum (above described).

The first swab has been withdrawn (sample number: 01).

The surface has been treated with vapour jet generated by Hygienio containing *ADANTIUM PLUS* disinfectant diluted at 2,5% for a period of about 60 seconds.

After a period of 5 minutes has been dragged the second swab (sample number: 07) in the right side of the contaminated surface.

8.2 Test execution mode

Microorganisms at 30°C

It was taken a sterile Petri plate and it was transferred through a sterile pipette 1 ml of the initial suspension (dilution 10⁻¹). Starting by the initial suspension, we have obtained some scalar dilutions until we arrived at the 10⁻⁴ dilution.

From every dilution has been withdrawn 1 ml, putting it in a Petri plate, where inside it was stilled for inclusion the PCA. Once solidified, the plates have been put to incubate at 30°C per 72 h ± 3 h. Then plates were been read through the use of colonies count, calculating the number of units forming colonies for ml.

8.3 Colonies count mode

The laboratory for the numeration techniques uses a plate for the dilution for almost two successive dilutions. After the period of necessary incubation, the colonies number is calculated.

MICROORGANISMS AT 30°C

These test methods include a maximum counting of colonies on the plate, equal to 300.

> For plates that don't contain colonies, the Laboratory expressed the result as <10 ufc per cm².

> For plates that contain a number of colonies between 1 and 3 on every sow plate, the Laboratory expresses the result as < 4 x d ufc per cm².

> For plates that contain a number of colonies between 4 and 300, the Laboratory applies, only after to have verified the proportionality of dilutions, the formula of the general case of ISO 7218, that is $N = \Sigma C / V \times 1,1 \times d$ where:

ΣC is the sum of the total number of colonies that were counted on every plate;

V is the volume of the inoculum stilled on every plate, in ml;

d is the dilution corresponding at the first considered dilution (d = 1 when referring to the liquid product not diluted of the sample test).

9. RESULTS

The results obtained are resumed in the following table and brought back in the certificate presented in the attachments (ALLEGATI VI- X)

MICROBIOLOGIC COUNT (UFC/ml)	
RESULTS	Pre Treatment
Total microbial charge control swab, control N01	$2,8 \times 10^3$ (UFC/cm ²)
	Post Treatment
Total microbial charge control swab, control N07 after a contact time of 5 minutes	< 1 (UFC/cm ²)

10. CONCLUSIONS

Results obtained with disinfection Hygienio- *ADANTIUM PLUS* system were very satisfactory, underlining the efficiency of the treatment that, during the contact time of 5 minutes, reconditions surfaces.

The test on microorganisms research turn out to be positive, satisfying completely the set goals.

As it can be noted from the test reports of the microorganisms count, the process reaches the highest efficiency in 5 minutes of contact, destroying completely the bacterial charge.

The achievement of this reduction in the time of a few minutes, let you continue quietly normal activities ambulatory.

This result is useful to the valuation of a economic return, in the field of Health operations, since it involves a reduction of the economic expense on materials, professional operators and mainly on the risk management.

11. ATTACHMENTS

- I. Bacteria test photos
- II. CERTIFICATION ACCREDIA LABORATORY JACARONI- "General requirements for the test and calibration laboratories competency".
- III. IQNet / DQS- Certification ISO 9001:2008 US LABS Srl.
- IV. INSTITUTE OF HEALTH-Approval of the system of warranty of production quality and/or sterilization of the US LABS Srl.
- V. TEST REPORT JACARONI DIAGNOSTIC 3656/2014.
- VI. TEST REPORT JACARONI DIAGNOSTIC 3662/2014.

ATTACHMENTS

Photo Bacterial Test

Foto 1: Bacterial Contamination of surface


Foto2 : Sampling swab before Hygienio treatment

Foto 3: Hygienio-Adantium Plus treatment

Foto 4: Sampling swab at 5 minuts treatment

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meets the requirements
of the standard

EN ISO/IEC 17025:2005 "General Requirements for the Competence of Testing
and Calibration Laboratories" standard

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The accreditation certifies the technical competence of the laboratory limited to the scope detailed in the attached Enclosure. The scope may vary in the time. The management system requirements in ISO/IEC 17025:2005 (Section 4) are written in a language relevant to Testing laboratories operations and meet the principles of ISO 9001:2008 and are aligned with its pertinent requirements.

The present certificate is valid only if associated to the annexed schedule, and can be suspended or withdrawn at any time in the event of non fulfillment as ascertained by ACCREDIA.

The in force status of the accreditation may be checked in the WEB site (www.accredia.it) or on direct request to appointed Department.

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
CERTIFICATO

IQNet

THE INTERNATIONAL CERTIFICATION NETWORK [®]

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US LABS S.r.l.

Strada di Recentino, 7
05100 Terni
Italy

has implemented and maintains a **Quality Management System**.

Scope:

Research, development, production and selling of products for disinfection and sterilization

Through an audit, documented in a report, it was verified that the management system fulfills the requirements of the following standard:

ISO 9001 : 2008

Date of certification 2011-12-02

Valid until 2014-12-01

Registration Number: 428803 QM08

Michael Drechsel
Managing Director of DQS GmbH
President of IQNet

IQNet Partners*:

AENOR Spain AFNOR Certification France AIB-Vinçotte International Belgium ANCE Mexico APCER Portugal CCC Cyprus
CISQ Italy CQC China CQM China CQS Czech Republic Cro Cert Croatia DQS Holding GmbH Germany DS Denmark
ELOT Greece FCAV Brazil FONDONORMA Venezuela ICONTEC Colombia IMNC Mexico
Inspecta Certification Finland IRAM Argentina JQA Japan KFQ Korea MSZT Hungary Nemko AS Norway NSAI Ireland
PCBC Poland Quality Austria Austria RR Russia SII Israel SIQ Slovenia SIRIM QAS International Malaysia
SQS Switzerland SRAC Romania TEST St Petersburg Russia TSE Turkey YUQS Serbia

IQNet is represented in the USA by: AFNOR Certification, CISQ, DQS Holding GmbH and NSAI Inc.

* The list of IQNet partners is valid at the time of issue of this certificate. Updated information is available under www.iqnet-certification.com

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APPROVAL OF QUALITY ASSURANCE SYSTEM FOR PRODUCTION AND/OR STERILIZATION

according to Annex V of EC Directive 93/42/EEC and following modifications
 (transposed in Italy by the DLgs n. 46 issued on 24.02.1997 and following modifications)

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*The Istituto Superiore di Sanità, Notified Body n° 0373, certifies that the quality assurance system for the
 production and/or sterilization enforced by*

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Is in compliance with the applicable requirements of EC Directive 93/42/EEC and following modifications.

Certificato n° Certificate no.	066	QPZ	1296	11	Addendum n° addendum no.	//	//
-----------------------------------	-----	-----	------	----	-----------------------------	----	----

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U.S. LABS S.r.l.