

Steam Sterilizers
Matachana S1000 RBE



Ideal solution for potentially infectious waste sterilization

The high versatility and operational simplicity of the S1000 RBE turns into the ideal sterilizer to treat healthcare waste or contaminated samples in a wide range of sterilizing centres, hospitals, microbiological laboratories, universities and research centres.

The MATACHANA Series S1000 RBE is a medium/high production capacity sterilizer that comprises several models: chamber volume between 445 and 893 liters, one or two doors, with or without steam generator.

Features

- Chamber and one-section-only jacket made of high-quality stainless steel AISI 1.4404 (AISI 316L).
- Automatic-pneumatically activated vertical sliding doors with safety Interlock and stop button.
- Control by industrial micro-computer and user friendly colour touch screen. Alphanumeric printer for cycle data recording.
- Stainless steel connection of steam pipes, 1.4404 (AISI 316L) quality.
- Vacuum by means of ejector (Venturi system): low maintenance and noiseless.
- Independent water tanks for vacuum system and steam generator, with cooling coil for recovering heat from condensates. Water consumption savings (when the external cooling system option is installed).
- Mechanized door gasket groove (longer gasket lifespan and easy maintenance).
- Air treatment system. It is extracted from the chamber through an absolute HEPA filter, with condensate treatment inside the chamber, to ensure a correct contamination-free sterilization.

Programs

- 2 test programs (Vacuum Test & Bowie & Dick Test),
- 2 programs specially designed for the treatment of medical waste at 134 °C: with 15 and 30 minutes sterilization, in accordance with configuration and load type.

Models	1006VR-1 1006ER-1	1006VR-2 1006ER-2	1008VR-1 1008ER-1	1008VR-2 1008ER-2	1010VR-1 1010ER-1	1010VR-2 1010ER-2	1012VR-1 1012ER-1	1012VR-2 1012ER-2
Technical specifications	1 door	2 doors						
Total dimensions (mm) Height, excluding filter (*) Width Depth	1954 996 1314	1954 996 1336	1954 996 1614	1954 996 1636	1954 996 2054	1954 996 2076	1954 996 2313	1954 996 2338
Chamber dimensions (mm) Height Width Depth	670 670 998	670 670 996	670 670 1265	670 670 1296	670 670 1735	670 670 1733	670 670 2000	670 670 1998
Power (kW) WITH generator (ER) WITHOUT generator (VR)	51 3	51 3	64 3	64 3	64 4	64 4	64 4	64 4
Chamber volume (liters)	445	445	565	578	774	773	893	892
Load capacity/cycle (kg) (**)	45-75	45-75	65-105	65-105	80-135	80-135	90-155	90-155



Aware of the problems which hospitals and research centers are currently facing, MATACHANA has developed a system which integrates itself within the special biocontaminated waste circuit.

- Classification and collection of the waste in suitable receptacles at the places where it is generated.
- Transportation flow within the hospital and its treatment in the sterilizer have all been thoroughly studied to elaborate a system that is simple and easy to install.
- The sterilization process converts the special waste into urban grade waste.
- The post sterilization process reduces the volume of waste treated.

All the power of technology in sterilization

- Working with the sterilizer is very easy. The user operates from the touch screen. Starting from the main menu, the screen guides the user, through the Submenus, toward the operations required.
- The follow-up of the cycle can be carried out directly from the screen, given that it allows a graphic visualization of the process in "real time".
- The control unit can integrate the Teleservice option ("on line" remote maintenance) and is provided with access doors to connect it to the PC and the bar code reader (optional).

Environmental Management 💙

Our constant eagerness for the technological innovation aimed at the sustainability and improvements of our products has taken us to control minutely the life cycle of our devices as well as their environmental impact. Such a control has allowed us to improve the manufacturing process of our devices obtaining considerable energetic saving improvements with regards to previous series.

••• matachana





I. I. 🗰 🗲

www.matachana.com



MATACHANAGROUP

Innovation is the way forward

Since our foundation 50 years ago, our mission has been to provide the best service, bringing our knowledge and field experience to our customers to facilitate their daily work, allowing them to be efficient in the production whereas keeping rigorously the quality.

Our history has taken us from the beginnings in 1962 to the transformation in 2000 in MATACHANA GROUP, with its subsidiaries in Germany, France, USA, Argentina, Malaysia and China, as well as the collaboration of direct distributors in over 110 countries.

In MATACHANA GROUP we are aware that Training and Service contribute decisively to achieve customer satisfaction. For this reason, we have always invested in the development of these two areas of activity, which enables a direct contact with customers and help us to develop together a continuous improvement process.

Technical Engineering Support

An assistance provided by engineers, highly skilled expert technicians and support staff, all committed to ensure the proper equipment operation and condition.

Training Academy



Committed to the responsibility to obtain the optimum operation from our equipment, in MATACHANA GROUP we invest in the training of future users.

In addition, we offer training courses regularly in the 5 continents, providing this service to all our customers.

Environmentally friendly \heartsuit

The sterilizers are designed and manufactured using the latest technologies on the market to achieve the best results in terms of energy savings and reduction in water consumption.

Our Production Center complies with the ISO 14001 Environmental Management System and ISO 50001 Energy Management.

Quality

The MATACHANA devices are developed, manufactured and tested within a strict quality control according to the International Standard ISO 9001.



MATACHANAGROUP Headquarters · Almogàvers,174 E-08018 Barcelona

Spain Tel. (+34) 93 486 87 00 hospitales@matachana.com International Tel. (+34) 93 300 85 04 int@matachana.com